


Inoculate replicate cultures,  
grow to mid-log phase,  
normalize cell densities

1



## Treatment

Treat with sublethal  
( $2 \times IC_{50}$ ) antibiotic  
for 16 hours







Plate on agar with  
sublethal ( $2 \times IC_{50}$ ) antibiotic,  
perform ScanLag, and  
cluster colony types



Pick colonies,  
phenotype, and  
identify mutations




Whole genome sequencing



2


Plate untreated controls  
on antibiotic-free agar



## Controls

Perform ScanLag  
on untreated controls

3



4